



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Plasma Fission (New)**

More information about Plasma Fission at [www.plasmafission.com](http://www.plasmafission.com)

- New,
- Made in Poland,
- A device for performing treatments using a plasma beam - it has many features that distinguish it from other devices,
- An excellent alternative to the scalpel. Skin lesions can be excised with a scalpel, but it is much better to treat them with a plasma beam and remove them in this way. Another purely aesthetic issue is wrinkle reduction and blepharoplasty. No other method has the same effect as using tissue coagulation with a plasma beam,
- The device is CE medical certified,
- Main applications:
  - Elimination of molluscum contagiosum,
  - Elimination of nevi (e.g. Miescher's),
  - Removal of seborrhoeic warts,
  - Removal of fibromas,
  - Removal of yellow tufts,
  - Alleviating hyperkeratosis,
  - Removal of hemangiomas,
  - Removal of telangiectasias,
  - Reduction of pigmented lesions,
  - Supporting dermatochalasis (blepharoplasty),
- Characteristics:
  - The plasma generated at the interface between the electrode of the device and the patient's skin produces a thermal effect and consequently coagulates the protein and vaporises the tissue,
  - The procedure does not cause bleeding, as the coagulated vessels are immediately closed,
  - It is almost painless and safe for the patient,
  - The high precision of the device allows the operator to control the field of action and the depth of penetration of the beam, making Plasma Fission a great alternative to the scalpel,
  - Plasma Fission is fully wireless,
  - No passive electrode,
  - Battery powered,
  - Clear battery level indicator,
  - Magnetic charging connector in the device protects it from breaking and facilitates quick connection of the charger,
  - Energy saving - when not in use, the unit will sleep and signal reduced power consumption and switch off after a period of time,
  - The aluminium carrying case will protect Plasma Fission during transport,
- Specifications:
  - Ergonomic shape to hold the device securely - weight is only 135 grams,
  - Battery powered,
  - Intuitive operation - clear OLED colour display,
  - Power control,
  - Plasma generator with automatic voltage regulation ensures constant beam power and safety during treatment,
  - Extensive software and hardware protection system,
  - MagClick™ electrode attachment system for quick electrode replacement,
  - Disposable sterile treatment electrodes with a diameter of 0.15mm or allow precise targeting of the plasma beam,
  - Built-in 8Wh rechargeable battery allows performing treatments all day without recharging,
- Valid technical passport issued,
- Warranty 24 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

