



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Clamp for medical rails with vertical lamp outlet (New)

- New,
- Made in France,
- Fits all medical rails,
- Ideal for use on a gynecological chair equipped with a handrail in the seat,
- Made of anodized aluminum,
- With mounting plate for mounting the lamp,
- Supplied with lamp adaptor plate. Tightening with allen wrench (included) + lock washer and nut,
- Maximum opening: 60 mm,
- Weight: 400 g,
- Warranty: 12 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

