



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Clamp for medical rails with horizontal lamp outlet (New)**

- New,
- Made in France,
- Fits all medical rails,
- With mounting plate for lamp mounting,
- Supplied with lamp adaptor plate. Tightening with allen wrench (included) + lock washer and nut,
- Made of anodized aluminum,
- Weight: 400 g,
- Warranty: 12 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

