



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Mammographic monitors NDS DOME E5 (Reconditioned)

Used, very good condition,

U.S.A. production,

Power supply: through a 12 V power supply,

Frequency: 50/60 Hz,

Power: 72 W,

The Dome E5 grayscale display is designed for the most demanding diagnostic applications, including mammography and radiology, both vertically and horizontally,

The Dome E5 display provides smooth image manipulation and window alignment, and supports open architecture graphics standards for 3D and volumetric rendering applications,

With full DICOM compatibility, automatic calibration and remote monitoring, the Dome E5 is ideal for the most demanding PACS systems,

Technology: TFT AMLCD Dual Domain IPS, 3 subpixels per pixel,

Screen size (diagonal): 21.3" (541 mm),

Display resolution (pixels): 2560 x 2048,

Pixel pitch: 165 μ m, 154 dpi,

Angle of vision: 170° horizontally and vertically,

Response time: 50 ms,

Brightness: 1000 cd / m² (default 500 cd / m²),

Contrast ratio: 800:1,

DICOM LUT: 3061 unique shades of grey,

Video input: Dual Link DVI-D,

Data sizes in pixels: 8-bit and 10-bit grey scale,

Compatible with VGA: 640 x 480 to 1280 x 1024,

USB interface: built-in HUB USB 2.0 on the rack,

The desk base includes height, horizontal / vertical rotation and tilt adjustment,

The set includes:

Two NDS Dome E5 monitors,

NVIDIA Quadro K2000D graphics card,

Display size without base: 470 x 387 x 103 mm,

The weight of a monitor without a base: 6.4 kg,

Monitor weight with base: 10 kg,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan

