



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Video - OLYMPUS gastroscopes GIF-H190 (Reconditioned)

Used, very good condition,
Japanese production,
Excellent image quality and maneuverability in a slender design,
HDTV clarity together with Narrow Band Imaging (NBI) provides better visualization and the water jet channel helps for clearer observation during procedures,
Field of vision: 140°,
Direction of vision: Forward view,
Deep down: 2 - 100 mm,
Outer diameter of the distal end: 9.2 mm,
Outer diameter of the inlet pipe: 9.2 mm,
Working length: 1030 mm,
Internal channel diameter: 2.8 mm,
Minimum visible distance: 3,0 mm from the distal end,
The bending angles of the tip:
Up: 210°,
Down: 90°,
Right: 100°,
To the left: 100°,
Total length: 1350 mm,
Compatibility with EVIS EXERA III systems (190),
Original suitcase included,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

