



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Video - OLYMPUS colonoscope CF-H190L (Reconditioned)**

Used, very good condition,  
Japanese production,  
Excellent HD clarity,  
A brighter narrowband imaging,  
The rapid response technology and water jet make the CF-H190L a robust tool for routine colonoscopy,  
HDTV and narrowband imaging enable optical diagnostics and provide access to DISCARD strategies,  
Reduced risk of water ingress with fully waterproof One-touch connector,  
Field of vision: 170°,  
Direction of vision: Forward view,  
Deep down: 2 - 100 mm,  
Outer diameter of the distal end: 13.2 mm,  
Outer diameter of the inlet pipe: 12.8 mm,  
Working length: 1680 mm,  
The inner diameter of the channel: 3.7 mm,  
Minimum visible distance: 3,0 mm from the distal end,  
The bending angles of the tip:  
Up: 180°,  
Down: 180°,  
Right: 160°,  
To the left: 160°,  
Overall length: 2005 mm,  
Compatibility with EVIS EXERA III systems (190),  
Original suitcase included,  
He has a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (Poland only): Instalments, Leasing, Loan

