



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MCM 404 Peristaltic pump (Reconditioned)

Used, very good condition,
German production,
Power supply: 230 V and rechargeable,
Frequency: 50/60 Hz,
Power consumption: 9 VA,
Scope of delivery: 1 - 999 ml/h,
Graduation: 1 ml,
Timer set to: 99 hours and 59 minutes,
Flow accuracy: +/- 5%,
Closure and volume alarm for the bolus,
Included is a drop (flow) detector on a spiral tube,
Dimensions: 16 x 22.5 x 22.5 cm,
Weight: 3.8 kg,
User manual in German (PDF),
Cat 02,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

