



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Peristaltic pump MCM 404 (Reconditioned)**

Used, very good condition,  
German production,  
Power supply: 230 V and rechargeable,  
Frequency: 50/60 Hz,  
Power consumption: 9 VA,  
Scope of delivery: 1 - 999 ml/h,  
Graduation: 1 ml,  
Timer set to: 99 hours and 59 minutes,  
Flow accuracy: +/- 5%,  
Closure and volume alarm for the bolus,  
Included is a drop (flow) detector on a spiral tube,  
Dimensions: 16 x 22.5 x 22.5 cm,  
Weight: 3.8 kg,  
User manual in German (PDF),  
Cat 03,  
He has a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (Poland only): Instalments, Leasing, Loan

