



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Orthopaedic support for the leg (Reconditioned)

Used, good technical condition,
German production,
Support made of soft foam covered with organic leather in white,
The upholstery has light tears,
Support ended in a stiffener,
The support ensures the lifting of the limb in the treatment of venous system diseases, e.g. inflammations and vein thrombosis,
The splint can also be used after surgeries and orthopaedic injuries for functional fracture treatment and post-traumatic rehabilitation,
The rail can be used by placing an angled platform under the thigh so that the knee and foot rest on the hill,
In addition, the support has a foot brace,
Easy to keep clean,
Dimensions:
Length: 77 cm,
Height: 46 cm,
Width: 24 cm,
Weight: 2 kg,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

