



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Diathermy VALLEYLAB FORCE 300 with argon (Reconditioned)

Used, very good condition,

U.S.A. production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 850 VA,

Foiled control panel,

Clear digital displays,

It has the following functions:

Cutting - 300 W,

Coagulation - 120 W,

Bipolar - 70 W,

Force GSU System argon diathermy included:

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 40 W,

Flow range:

Laparoscopic mode: 0.5 - 4 l/min,

Open mode: 2 - 16 l/min,

Acoustic and visual low pressure alarm,

Gas flow control,

Argon diathermy placed on a table with space for two argon bottles (bottles included) and two shelves for accessories with dimensions: 44 x 26 cm,

In photo number 10 you can see a wearable disposable neutral and monopolar electrode,

Dimensions of diathermy: 40 x 46 x 14 cm,

The dimensions of the set: 46 x 78.5 x 100 cm,

The weight of diathermy: 7 kg,

The weight of the set: 63 kg,

Volume: 65 dB,

The operating instructions are in German,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan

