



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Liposuction compressor EKOM DK50 2V/M (without box, with membrane dryer) (New)

A new one,
Slovakian production,
Technical data:

Rated voltage: 230V/ 50 (60) Hz

Oil-free,

Frequency: 3x400V/ 50 Hz

Power: 1.1 kW, 1.2 kW (3x400/50)

For power: two units,

Number of dental sets: 2

Compressor capacity (working pressure 5-7 bar):

108 l.min- 1/ 5 bar

Tank capacity: 25 litres.

Operating mode: S1-100%

The volume:

≤ 69 dB [A]

Weight: 63 kg/ 61 kg (3x400/50),

Degree of air drying: dew point -20°C

For power: two units,

Dimensions (Sz x g x w): 460x565x720 mm,

He has a current Technical Passport issued,

Warranty 12 months,

Possibility of financing (Poland only): Instalments, Leasing, Loan

