



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Filing cabinet - FEBRU (Reconditioned)

Used, very good condition,
German production,
Cabinet made of MDF boards,
Drawers made of powder-coated steel,
Drawer front panels made of plastic,
The cabinet is equipped with 4 drawers:
The first type of tray with depth dividers: 2 cm,
Three more drawers with dimensions: 33 x 70 x 19 cm,
Single-row full extension drawers,
Only one drawer can be opened at a time,
Comfortable handle for opening drawers,
Dimensions: 43 x 80.5 x 72.5 cm,
Weight: 42 kg,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

