



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Filing cabinet - FEBRU (Reconditioned)**

Used, very good condition,  
German production,  
Cabinet made of MDF boards,  
Drawers made of powder-coated steel,  
Drawer front panels made of plastic,  
The cabinet is equipped with 4 drawers:  
The first type of tray with depth dividers: 2 cm,  
Three more drawers with dimensions: 33 x 70 x 19 cm,  
Single-row full extension drawers,  
Only one drawer can be opened at a time,  
Comfortable handle for opening drawers,  
Dimensions: 43 x 80.5 x 72.5 cm,  
Weight: 42 kg,  
He has a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (Poland only): Instalments, Leasing, Loan

