



10 L oxygen therapy kit with trolley (Reconditioned)

Cylinder capacity: 10 litres,
Working pressure PW: 200 bar,
Test pressure PH: 300 bar,
Equipped with a 3/4" valve,
The cylinder is filled with medical oxygen: 1,6 m3,
It is certified and validly legalized,
Included in the set is an oxygen regulator with flow control and humidifier (type of regulator available depending on stock) and a choice:
Adult oxygen mask with a 2.1 m draught tube,
Pediatric oxygen mask with a 2.1m drain,
Adult nasal cannula with a 2.1 m drain,
The set includes a trolley for transporting cylinders:
Stable and durable,
Designed for the carriage of cylinders with max. diameter. 15 cm,
Made of powder-coated steel,
Rubber, non-slip mat on base,
Big, comfortable handle,
Two wheels with a diameter: 12.5 cm,
Cylinder dimensions:
Height: 91 cm,
Diameter: 14 cm,
Pictures,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

