



## **10 L oxygen therapy kit with trolley (Reconditioned)**

Cylinder capacity: 10 litres,  
Working pressure PW: 200 bar,  
Test pressure PH: 300 bar,  
Equipped with a 3/4" valve,  
The cylinder is filled with medical oxygen: 1,6 m3,  
It is certified and validly legalized,  
Included in the set is an oxygen regulator with flow control and humidifier (type of regulator available depending on stock) and a choice:  
Adult oxygen mask with a 2.1 m draught tube,  
Pediatric oxygen mask with a 2.1m drain,  
Adult nasal cannula with a 2.1 m drain,  
The set includes a trolley for transporting cylinders:  
Stable and durable,  
Designed for the carriage of cylinders with max. diameter. 15 cm,  
Made of powder-coated steel,  
Rubber, non-slip mat on base,  
Big, comfortable handle,  
Two wheels with a diameter: 12.5 cm,  
Cylinder dimensions:  
Height: 91 cm,  
Diameter: 14 cm,  
Pictures,  
He has a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (Poland only): Instalments, Leasing, Loan

