



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Operating table FAMED LIVE SU-06.0 (Reconditioned)

Used, very good condition,
Upholstery's good,
Polish production,
Power supply: 230 V or rechargeable,
Frequency: 50/60 Hz,
Power: 300 W,
8-segment,
Designed for operations and procedures in general surgery, vascular surgery, cardiac surgery, neurosurgery, urology, gynaecology, proctology, laryngology, ophthalmology, endoscopy, laparoscopy, trauma surgery, oncology, dental, plastic surgery and others,
Made of stainless, acid resistant steel,
The chassis has a central wheel lock,
The top has removable, seamless, antistatic mattresses resistant to disinfectants,
Electrical adjustment with the multifunctional remote control:
Heights in range: 79 - 117 cm,
Backrests in the range -24° to +78°,
Foot segment in the range from -45° to +80° (single leg adjustment possible),
Trendelenburg's position in the range: 26°,
Anti-Trendelenburg's position in the range: 28°,
Lateral tilt in range: 18°,
Longitudinal thrust in the range: 28 cm,
Central brake,
The "0" button,
Manual adjustment by means of gas cylinders:
Headrest segment,
The lumbar section,
The shank segment,
Leg segments are sideways reversed,
The tabletop segments are permeable for X-rays and allow X-rays to be taken and the patient to be monitored with the C-arm,
The set includes a trolley for transporting FAMED ŻYWIEC WB-06.1 tops with central brake and directional lock,
Operating table dimensions:
Sleeper width: 52 cm,
Length of the bed: 213 cm,
Thickness of mattress: 7 cm,
Total width: 64 cm,
Total length: 213 cm,
Total height: 79/117 cm,
Wheel diameter: 15 cm,
Dimensions of the transport trolley:
Width: 86.5 cm,
Length: 145 cm,
Height: 59/79 cm,
Wheel diameter: 20 cm,
The weight of the operating table: 273 kg,
The weight of the transport trolley: 70 kg,
Version A,
Cat 02,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

