



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Respirator - ventilator assistant HEINEN + LOWENSTEIN SOMNIA 3 C-FLEX+ (Reconditioned)

Used, very good condition,
German production,
Power supply via 12 V power supply,
Frequency: 50/60 Hz,
Automatic adjustment of the therapeutic pressure to the patient's needs,
Compact design and very quiet operation,
The Smart RAMP function is adjustable from 0 to 45 minutes,
Automatic switch-on and switch-off of the device - the device detects whether the patient is wearing a mask and switches on or off automatically when the mask is removed,
The A-FLEX function reduces pressure during exhalation and at the same time mitigates pressure variations, minimizing effort during exhalation,
The SD card stores therapy data, allows you to assess the effectiveness of the therapy - no card included,
Internal memory that protects therapeutic data in case of damage to the memory card,
Information on AHI, CSR, snoring and leakage indices in the system,
Proactive auto-algorithm actively prevents the occurrence of (among others) apnea, shallow breathing and snoring,
Power supply unit with power cord included,
All placed in the original carrying bag,
The dimensions of the camera: 15 x 17 x 11 cm,
The weight of the set: 1.8 kg,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

