



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Progetti PG M9000 12.1 inch patient monitor (New)

A new one,
Italian production,
With a 12" LCD display, the monitor becomes the ideal device to take on the challenge of transporting and monitoring patients in critical condition,
This monitor can be connected to a monitoring station,
Ideal for any hospital ward
Equipped with the following modules:
ECG - 5 output ECG,
HR - heart rate measurement in the range 10~300 bpm
RESP - possibility to measure RESP breathing with chest impedance
NIBP - Measurement with an automatic oscillometer,
Digital SPO2 - measurement of SPO2,
PR - pulse measurement from 20 bpm to 250 bpm,
TEMP - measurement between 0.0°C - 50.0°C,
IBP - invasive blood pressure measurement, 2-channel, (can be retrofitted at additional cost)
Co2: (possibility of retrofitting for an additional fee)
Sidestream,
Mainstream,
Microstream,
Technical data:
Dimensions: 318x264x152mm,
Weight: ≤4.5 kg,
The power supply:
Power supply voltage: AC 100-240V, 50/60Hz,
Compliant with IEC 60601-1 and IEC 60601-1-2,
Recorder (printer for an additional fee):
Uses the thermal method for recording
Horizontal resolution: 16 points / mm (at 25 mm / s paper speed),
Vertical resolution: 8 points / mm,
Paper width: 50 mm
Paper length: 15 m,
Paper writing speed: 12.5/25/50 (mm/sec)
Maximum 3 tracks
Display:
Diagonal of the monitor: 12.1",
Display matrix: TFT LCD,
Resolution: 800x600 pixels,
A touch screen,
Battery
Type: Lithium rechargeable battery 11.1V/4.0Ah,
Working time under normal use and full charge: ≥ 240 minutes (2 batteries per 480 minutes),
Storage of data:
Trend data: Long trend: 168 hours, minimum resolution is 1 min High resolution trends: 2 hours, minimum resolution is 5 seconds ,
Parameter alarm event: 128 groups of parameter alarm events,
ARR event: 128 groups of ARR events and the associated course for each waveform. NIBP measurement result: 1000 groups,
Included:
Monitor,
Reusable SpO2 sensor for adults,
Adult cuff (25 cm, 35 cm),
NIBP tube (3 m) with connector,
5- drain snap-in connector ECG cable,
A sealed acid battery,
Power cable,
Ground wire,
The manual,
CE marking according to MDD93/42/EEC,
He has a current Technical Passport issued,
Warranty - 12 months,
Possibility of financing: Installments, Leasing, Loan.

