



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Signal monitor PROGETTI PG S10 (New)

A new one,
Italian production,
The monitor is equipped with an 8-inch TFT colour LCD screen (the touch screen is optional for an additional fee),
Portable lightweight and robust design makes it easy to move the monitor,
Thanks to the flexible configuration of parameters it is used in various medical wards,
The lithium-ion battery provides up to 12 hours uninterrupted operation,
The large and coloured font display setting makes it easy to monitor the patient,
Selectable mode for adult, pediatric and neonatal patients Wired / wireless CSM, HL7 to HIS support, optional graphical and tabular trend overview
48h full waveform overview for each patient (stored on SD cards)
Features and benefits of the PG S10 monitor,
8-inch TFT LCD colour screen (touch screen at extra charge),
Portable, lighter and robust design,
Flexible parameter configurations for different critical environments,
Lithium-ion battery (up to 12 hours uninterrupted operation),
Set to display large and colorful font,
Spot check and continuous monitoring mode,
Choice for adults, children and newborns,
Wired / wireless CSM, HL7 to HIS protocol support (optional extra charge),
Graphical and tabular overview of trends,
48-hour full review of test results for each patient (stored on SD cards),
Thermal printer (optional extra charge),
Equipped with the following modules:
NIBP - Measurement with automatic oscillometers,
SPO2 - digital measurement of SPO2,
TEMP - Measuring range: 5-50 ° C (41 ~ 122 ° F),
PR - pulse measurement from 30 bpm to 254 bpm,
Technical data:
Power supply: AC 100-250 V, 50/60 Hz,
The working environment:
Temperature: 5-40°C
Humidity: <85%
Patient range: adults, children, newborns
You can choose men in English, Spanish, Russian, French, Turkish, German, Polish, Romanian,
Included:
Monitor
SpO2 sensor for adults
NIBP cuff for adults
Skin temperature probe
He has a current Technical Passport issued,
Warranty: 12 months,
Possibility of financing: Installments, Leasing, Loan.

