



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PROGETTI EPG 6 View Plus unit (New)

A new one,
Italian production,
The PROGETTI EPG 6 View Plus is the classic PROGETTI electrocardiograph, an icon of reliability and precision,
EPG 6 VIEW PLUS is a 6/12 channel device and is ideal for the general practitioner, but also for the specialist and his or her needs,
Thanks to the compact keyboard, this device will allow the doctor to create a database in real time and transfer the data to his computer,
Parameter calculation and HES interpretation make EPG 6 VIEW PLUS the ideal electrocardiograph in any situation,
Features and benefits:
3/6/12 channel electrocardiograph with analysis and interpretation,
Internal memory for over 300 tests,
Pacemaker pulse detection (pacemaker),
Display: 5.7 " TFT LCD colour display,
USB output,
Computer software for data management,
A one-time record: 3/6/12 channels,
Record width: 112 mm,
External USB printer: A4 printout (210x297),
Alphanumeric keypad,
Automatic isoelectric line adjustment,
Constantly measuring heart rate (HR),
Acoustic signaling of detected stimuli,
Easy loading of paper into the printer,
Save the test copy to PENDRIVE, in EN1064 (ECG-SCP),
Signalling incorrect connection to individual electrodes,
Technical data:
Display: 5.7" graphics, colour TFT 320x240,
Power supply: 7.2V lithium-ion battery; 2200 mAh Li-ion,
Power supply: 90-240 V, 50/60 Hz,
Power consumption: <30VA,
The frequency band: 0.05: 150 Hz,
CMRR:> 100 dB ,
Entry impedance: $\geq 10 \text{ M}\Omega$,
Number of drains: 12 standard or Cabrera,
Sensitivity: 2.5 mm / mV; 5 mm / mV; 10 mm / mV; 20 mm / mV $\pm 5\%$
Recording speed: 5; 10; 25; 50 mm/sec,
Digital interference filtration: 25.35.50.60 Hz
Analysis and interpretation: HES according to EN 60601-2-51,
Sampling rate: 1000 Hz,
Converter: 12bit,
Resolution: 0.25uV / bit,
Dynamic range: 10 mVpp,
Dimensions: 260 W x 52 H x 220 D mm,
Weight: 1.8 kg,
Operating temperature: + 10°C ~ 40°C,
Operating humidity: 25% ~ 95% non-condensing,
Standards of compliance: EN 60601-1, EN 60601-1-2, EN 60601-2-25, EN 60601-2-51,
Safety: Type of protection CF (EN60601-1) Class I Class / Group: Class A / Group 1 (CISPR-11),
Included:
EKG device,
A 10-guided patient cable,
A set of electrodes,
A roll of paper,
Power cable,
Gel flacon ECG 260 ml,
A quick guide,
He has a current Technical Passport issued,
Warranty: 12 months,
Possibility of financing: Installments, Leasing, Loan.

