



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **PROGETTI EPG 12 View unit (New)**

A new one,  
Italian production,  
High-end electrocardiograph, thanks to the 8.4 inch touch screen,  
This device meets the expectations of every cardiologist,  
Main features: data storage, report management software, interpretation and various visualisations,  
EPG 12 VIEW is designed for those who want electrocardiography to be a distinctive feature of their work,  
Features and benefits:  
Display: 8.4 inch TFT touch screen,  
Preview before printing,  
Simultaneous recording of 12 ECGs,  
Automatic measurement and interpretation,  
120 seconds 12 discharges, recorded mileage,  
Up to 120 types of arrhythmia analysis,  
Large buttons, alphanumeric keypad,  
Built-in long-life lithium battery,  
Different print formats, support A4 210 mm, roll 216 mm,  
USB output (at extra charge),  
Technical data:  
Display: colour LCD 800 x 600 pixels,  
Size: LCD: 8.4 inches,  
Power supply: built-in lithium battery, 14.4V,  
Power supply: AC: 100 V-240 V, 50/60 Hz,  
Configuration: 12 channels  
The frequency band: 0.05: 150 Hz,  
CMRR: > 105 dB ,  
Entry impedance:  $\geq 50 \text{ M}\Omega$ ,  
Number of drains: 12 standard,  
Defibrillation protection: internal,  
Time constant:  $\geq 5 \text{ s}$   
Calibration voltage: 1 mV,  $\pm 1\%$   
Sensitivity: 2.5 mm / mV; 5 mm / mV; 10 mm / mV; 20 mm / mV ,  
Noise level:  
AC Filter: 50/60 Hz - closed,  
EMC filter: 25/35/45 Hz / closed,  
Drift filter: 0.05/0.15/0.25/0.50 Hz,  
Low pass filter: 70/100/150 Hz,  
Recording mode: Auto / Manual,  
Record format: 12x1, 3x4, 6x2, 3x4 + 1R, 3x4 + 3R, 6x2 + 1R,  
Automatic measurement function and automatic diagnosis,  
It can analyze up to 122 types of arrhythmias,  
Recorder: thermal printer,  
Roll of writing paper,  
Recording speed: 5; 10; 25; 50 mm/sec,  
Input:  $\geq 100 \Omega$ : Sensitivity 10 mm / V  $\pm 5\%$ ,  
Output:  $\leq 100 \Omega$  Sensitivity 1 V / V  $\pm 5\%$ ,  
Memory: Standard 300 ECG,  
ECG preview,  
Connect to an external printer (only in automatic mode),  
Dimensions: 410x316x136 mm,  
Included:  
EKG device,  
A 10-guided patient cable,  
A set of electrodes,  
A roll of paper,  
Charger,  
Gel flacon ECG 260 ml,  
Operating instructions  
He has a current Technical Passport issued,  
Warranty: 12 months,  
Possibility of financing: Installments, Leasing, Loan.

