



DRAGER OXYLOG 2000 ambulance respirator - ventilator (Reconditioned)

Used, very good condition,
German production,
Power supply via 12 V power supply or 1000 mAh battery,
Built-in manometer: -10 to 80 mbar,
Liquid crystal display,
Language of the device: German,
Adjustable ventilator functions:
Ventilation frequency (Freq): 0 - 40 1/min,
Tidal volume (Vt): 0.1 - 1.5 L,
The ratio of the time of inhalation and exhalation (Ti:Te): 1:3 / 1:2 / 1:1,5 / 1:1 / 2:1,
Maximum airway pressure (Pmax): 20 - 80 mbar,
End-of-expiratory pressure (PEEP): 0 - 15 mbar,
Ventilation mode switch: SIMV/CPAP or IPPV/SIPPV,
Mixing switch: 100% O2 or 60% O2,
The set includes:
A set of respiratory tubes,
Respiratory valve (valve), flow sensor, angular connector,
Oxygen supply hose with AGA tip,
Patient mask size 5,
Power supply with power cord,
User manual in Polish (PDF),
Dimensions: 21.5 x 21 x 17 cm,
The weight of the set: 6.1 kg,
Cat O2,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

