



### **Hospital bed INVACARE 750.420.35 (Reconditioned)**

Used, very good condition,  
American production,  
Power supply 230 V,  
Frequency 50-60 Hz,  
Power: 70VA,  
Very quiet engines,  
Electrically adjustable in many planes: up-down, headrest, footrest, sitting position, thigh backrest segment adjustment,  
The bed has four-segment beds made of aluminium,  
All bed functions are electrically controlled with the remote control,  
Metal bed frame, powder-coated (resistant to bruising),  
Gables and sides made of furniture board,  
Side railings folding along the bed frame,  
Maximum load 220 kg,  
Mattress size: 198 x 82 x 12 cm,  
General dimensions: 213 x 101 x 69/115 cm,  
Mobile (wheels with central brake),  
He has a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
3 month for the international market,  
Possibility of financing (Poland only): Instalments, Leasing, Loan

