



Hospital bed INVACARE 750.420.35 (Reconditioned)

Used, very good condition,
American production,
Power supply 230 V,
Frequency 50-60 Hz,
Power: 70VA,
Very quiet engines,
Electrically adjustable in many planes: up-down, headrest, footrest, sitting position, thigh backrest segment adjustment,
The bed has four-segment beds made of aluminium,
All bed functions are electrically controlled with the remote control,
Metal bed frame, powder-coated (resistant to bruising),
Gables and sides made of furniture board,
Side railings folding along the bed frame,
Maximum load 220 kg,
Mattress size: 198 x 82 x 12 cm,
General dimensions: 213 x 101 x 69/115 cm,
Mobile (wheels with central brake),
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
3 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

