



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Hospital bed PRAISTECH 1 (New)

New,
Polish production,
The bed has a stable steel construction with a fixed height coated with powder paint,
Two-segment beds are filled with steel mesh also powder coated,
The tops of the bed are filled with damage resistant and easy to clean HPL,
The structure has 4 wheels with brakes,
For a surcharge, options such as
Mattress
Side folding barriers
Drip hanger
Shelf for patient's personal belongings
Hanger for patient information card
Technical data:
Total length: 2150 mm,
Total width: 980 mm,
Height (without mattress): 500 mm,
Sleeping width (with mattress): 850 mm,
Length of bed (with mattress): 2000 mm,
Recommended mattress dimensions: 2000x850x100 mm,
Elevation angle of the square segment: 0° to 70°,
Self-preservation of the back segment: By 110 mm,
Maximum load: 250 kg,
Wheel diameter: 125 mm,
He has a current Technical Passport issued,
12-month warranty.
Possibility of financing (Poland only): Instalments, Leasing, Loan

