



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CARL ZEISS MediLive ImageBox digital video recorder (Reconditioned)

Used, very good condition,

German production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power consumption: max. 50 W,

It is a digital device that captures images obtained by video cameras and records them on hard disk and archives them by "burning" them on CD-R disks in JPEG or TIFF formats of choice,

The device is very simple to operate, it is controlled by buttons on the control panel, an additional foot pedal, function buttons on the microscope hand griffins or a 14-function foot pedal of the operating microscope,

The MediLive ImageBox connects directly to the video source via standard Y/C (S-video) or BNC connectors,

Small size and weight of the device and simple operation make the ImageBox extremely useful for use in the operating room, clinic or doctor's office,

You can view the recorded images by attaching a monitor,

The source of the video signal can be e.g:

Cameras mounted on operating microscopes,

Cameras mounted on endoscopes,

Cameras mounted on slit lamps,

Other cameras used e.g. for home use,

The set includes:

MediLive ImageBox recording device,

Power cord,

The user manual in English,

Dimensions: 31.7 x 34.5 x 7 cm,

Weight: 4.5 kg,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan

