



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Infusion (syringe) pump FRESNIUS INJECTOMAT 2000 P (Reconditioned)

Used, very good condition,
French production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 23 VA,
Clear electronic graphic display and additional liquid crystal display of line pressure,
Choice of syringes: 20 and 50/60 ml,
Supply rate: 0.1 - 200 ml/h,
Bolus: 0.1 - 800 ml/h,
The infusion volume: 0.1 - 999.9 ml,
Acoustic alarms,
A syringe sized reader,
Programmable name of the drug,
Automatic bolus reduction after occlusion,
Display and setting of infusion pressure,
Quick start,
RS232 socket,
Comfortable carrying handle,
Dimensions: 32.5 x 19 x 11.5 cm,
Weight: approximately 2.2 kg,
User manual in German (PDF),
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

