



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Video - gastroscope OLYMPUS GIF-300 (Reconditioned)

Used, very good condition,
Japanese production,
Direction of vision: Forward view,
External diameter: 12.9 mm,
Working channel diameter: 3.2 mm,
Total length: 3000 mm,
Maximum bending angle of the endoscope tip up and down: 180° / 180°,
Max. Bending angle of the tip right and left: 160° / 160°,
Maximum angle of field of vision: 140°,
Fully compatible with processors from CV/CLV 100 to 190,
Waiting time of 14 to 60 days,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

