



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS EVIS EXERA II videogastroscope GIF-Q180 (Reconditioned)

Used, perfect condition,
Japanese production,
Direction of vision: Forward view,
External diameter: 8.8 mm,
Outer diameter of the tip: 8.8 mm,
Working channel diameter: 2.8 mm,
Working length: 1030 mm,
Total length: 1345 mm,
Deep down: 3 - 100 mm,
Maximum bending angle of the endoscope tip up and down: 210° / 90°,
Max. Bending angle of the tip right and left: 100° / 100°,
Maximum angle of field of vision: 140°,
Fully compatible with CV-160/140 processors,
Included is an original suitcase with dimensions: 59.5 x 43.5 x 13 cm,
Weight: 4.3 kg,
User manual in English (PDF),
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

