



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS Evis Exera II GIF-Q160 videogastroscope (Reconditioned)

Used, perfect condition,
Japanese production,
Direction of vision: Forward view,
External diameter: 9.2 mm,
Working channel diameter: 2.8 mm,
Working length: 1030 mm,
Total length: 1345 mm,
Deep down: 3-100 mm,
Max. Angle of the endoscope tip bend up and down: 210°/90°
Max. Bending angle of the tip right and left: 100°/100°,
Maximum angle of field of vision: 140°,
Fully compatible with CV-180/160/140 processors,
It's complete:
Original suitcase with dimensions: 60 x 43 x 13.5 cm,
Valves,
Lid,
Cat 01,
It has a current overview, and a Technical Passport,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan
OVERVIEW PICTURES OF THE DEVICE

