



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS CYF-V cystofiberoscope (Reconditioned)

Used, very good condition,
Japanese production,
Direction of vision: Forward view,
External diameter: 5.4 cm,
Working channel diameter: 2.2 mm,
Working length: 380 mm,
Maximum bending angle of the endoscope tip upwards: 210°,
Maximum downward bending angle of the endoscope tip: 120°,
Maximum angle of field of vision: 120°,
Fully processor-compatible: Olympus OTV-S7 (PAL) / Olympus OTV-SI (PAL) / Olympus EVIS EXERA II (PAL) / Olympus VISERA / VISERA Pro (PAL),
Included:
Cover for ETO gas sterilization,
Biopsy valves MAJ-579,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan,
OVERVIEW PICTURES OF THE DEVICE

