



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Infusion (syringe) pump FRESNIUS VIAL PROGRAM 2 (Reconditioned)

Used, very good technical condition,
French production,
Power supply: 230 V and via a built-in battery,
Frequency: 50/60 Hz,
Infusion pump for 2 syringes,
Possibility of supply from two syringes simultaneously,
Readable electronic graphic displays,
Flow rate ranges: 0.1 to 400 ml / h (default 200 ml / h),
Increase: 0.1 ml / h - adjustable flow rate range,
Automatic recognition of the syringe volume,
Flow rate: 50 to 500 ml / h with volume display,
Syringe selection: 10 - 50 ml,
Acoustic alarms,
RS232 socket,
Dimensions: 30 x 23 x 15 cm,
Weight: 4.4 kg,
Has a valid Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (only Poland): Installments, Leasing, Loan,

