



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



TAUTMANN space-saving wheelchair (Ex-Display)

A new one, post-exposure,
German production,
60% space saving by grouping the seats,
Stainless steel with special surface treatment,
Elegant look thanks to chromium-nickel coatings,
Foldable armrests,
Foldable footrest,
Ergonomic seat and backrest for maximum comfort,
Easy and comfortable to use thanks to the 360° swivelable wheels,
Central brake on the two rear wheels,
Dimensions:
Width: 690 mm,
Length: 1070 mm,
Height: 985 mm,
The width of the seat: 490 mm,
The depth of the seat: 455 mm,
The height of the backrest: 440 mm,
Backrest width: 270 mm,
Permissible load: 140 kg,
The weight of the chair: 18.5 kg,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

