



Docking station for up to 4 pumps B.BRAUN 8713140 (without handle) (Reconditioned)

Used, very good condition,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 80 VA,
The station has guides for infusion lines,
It has a place for mounting to a bracket - no bracket included,
Dimensions: 27.5 x 16 x 33 cm,
Weight: 3.3 kg,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

