



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Biochemical analyser SINNOWA B200 (Reconditioned)

Used, very good condition,
P.R.C. production. (11.2018r.),
Power supply: 230 V,
Frequency: 50/60 zHzHyHz,
Power: 550 VA,

The B200 automatically performs all the steps with the pipetting of the serum and reagents. Programming takes about 10 seconds. The instrument performs the recommended tests without user intervention, which saves time,

Full automatic - very fast and easy programming,

It has a system for cooling the reagents, which prevents their degradation,

Very high repeatability and precision of measurements,

The analyzer performs all biochemical tests,

Capacity: 140 tests per hour,

Very low cost per mark,

It has 64 items for samples and 20 for reagents,

Possibility to enter patient's data: gender, age, name and other notes,

It is possible to enter the office logo,

Memory of about a million results,

CITO samples can be inserted,

Possibility to preview the operation of the camera,

The set includes:

Biochemical analyzer B200,

Laptop with software,

Automatic puppet,

Reagent kit,

Full set of container trays,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan

