



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## ISO9001 Certified COVID-19 IgM/IgG Ab - Core tests (cassette - plate) (25 pcs) (New)

New,  
P.R.C. production,  
Core test is a certified cassette test for Covid - 19,  
Detects in human serum, plasma or whole blood the presence and type of antibodies (IgM and IgG) produced by the immune system of a human infected with the SARS - CoV - 2 coronavirus,  
Intended for professional in vitro diagnostics,  
The performance of the test has been tested and certified: medical CE, ISO 9001 and ISO 13485,  
IgM antibodies can be detected with the test after 4 - 12 days from infection, then IgG antibodies are detected,  
Approved by the Polish Office for Registration of Medicinal Products (URPL),  
Based on a detection method using colloidal gold, which allows the determination of IgM and IgG antibodies in patients,  
Reading of the result possible after 15 minutes,  
Specification:  
Specificity: 100 %, false positive 0 %,  
Sensitivity: 97%, false negative 2.4%,  
Product storage temperature: 2 - 30° C,  
Included:  
COVID test cassette - 19 IgM / IgG Ab - 25 pieces,  
Ampoule with buffer - 25 pcs,  
Pipette - 25 pieces,  
Instruction for use in Polish - 1 piece,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE

