



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **MEDELA fat transfer kit (Fat Collector / Fat Transfer) (New)**

New,  
Swiss production,  
The set includes:  
1000 ml fat container with a lid and overflow valve,  
Reusable stand made of stainless steel,  
Reusable silicone hose,  
Additional equipment for an extra charge:  
Disposable 17 cm hose with clamp and Luer Lock faucet,  
3000 ml fat container  
The container is made of high quality, unbreakable, impact-resistant polysulphone. It has a clear marking of liquid capacity. It contains no latex and is suitable for autoclave sterilization,  
The cover is made of reinforced polyamide (Santoprene gasket) with a conical patient connector and mechanical overflow protection. It does not contain latex and is suitable for autoclave sterilization,  
He has a current Technical Passport issued,  
Warranty 12 months,  
Possibility of financing (Poland only): Installments, Leasing, Loan.

