



### **Docking station for 6 pumps FRESENIUS KABI ORCHESTRA BASE A (Reconditioned)**

Used, very good condition,  
French production,  
Power supply: 230 V or via the built-in 6 V / 1.3 Ah battery,  
Frequency: 50/60 Hz,  
Power: 125 VA,  
The Base A docking station is a data hub and power supply for up to 6 Visio DPS modules and/or Visio DPS modules and/or MVP modules,  
It allows for faster therapeutic control,  
Communication interface compatible with the main PDMS software,  
Built-in RS232 interface,  
One single power cable for all connected modules,  
Battery life - minimum 2 hours,  
Equipped with a tripod with a comfortable carrying handle,  
No pumps included - it is possible to retrofit the FRESENIUS ORCHESTRA MODULE DPS+ infusion (syringe) pumps for an additional fee,  
Dimensions of the docking station:  
Height: 55 mm,  
Width: 315 mm,  
Depth: 200 mm,  
Height with tripod: 1015 mm,  
The weight of the station with the tripod: 4 kg,  
The weight of the docking station with pumps: 24 kg,  
He has a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (Poland only): Instalments, Leasing, Loan

