



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ALARIS GATEWAY workstation (Reconditioned)

Used, very good condition,

Swiss production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 460 VA,

The Alaris Gateway workstation has space for 5 ALARIS syringe pumps and 2 volumetric (volumetric) infusion pumps,

Today's intelligent infusion systems require intelligent communication and integration solutions. Carefusion products can integrate infusion data with the hospital's central clinical information system. These wired / wireless solutions provide total control and valuable infusion information 24 hours a day,

Numerous control and control parameters ensure that the Alaris Gateway Workstation provides detailed and precise information. Reliability is ensured by proven integration protocols,

The Alaris Gateway Workstation uses established standards to ensure high efficiency and performance to work reliably with your current system. This solution uses TCP/IP as a network protocol and HTTP as a protocol for basic communication,

The station is equipped with two height-adjustable drip-hangers,

Elimination of power cable bundles by using a single power cable,

Two mounting brackets built in,

No pumps included - it is possible to retrofit all ALARIS pumps from our offer for an additional fee,

Tripods on wheels visible in some pictures are not included in the set,

Dimensions:

Width: 38 cm,

Height: 101/149 cm,

Depth: 25.5 cm,

Weight: 12.6 kg,

User manual in English (PDF),

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan

