



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **NOUVAG Vacuson 40 suction pump / mammal (New)**

New,

The Swiss production,

The advantages of the VACUSON 40 Mammal 40:

Economical thanks to excellent price/performance ratio,

Precise and infinitely adjustable suction capacity,

Robust, reliable and proven technology,

Low noise and vibration levels,

Oil-free piston pump,

Vacuson 40 devices are designed for suction and removal of fluids and intraoperative contents with adjustable suction force, which allows the device to adjust its operation to the needs of the treatment,

The suction cup has a manometer to control the settings,

Piston technology used in Vacuson 40 mammals allows to obtain high vacuum parameters with pressure up to 0.9 bar and to work for many years without any decrease in parameters,

Thanks to their advantages, the devices can be used in both small and very large procedures,

The Vacuson 40 is equipped with an automatic overfilling system,

Standard equipment:

Control unit VACUSON 40 (pressure gauge 0-1 bar, vacuum regulator, 2 jar holders, 3 sockets: for ON/OFF footswitch, foot regulator and power cord, suction pipe with set of connection wires and disposable bacterial filter),

Footswitch On/Off,

Disposable bacterial filter (10pcs),

Power cable (5m),

Additional equipment (surcharge):

Jars of different capacities,

The "STAND" mammal trolley (2 shelves)

SANICAR mammal trolley (3 shelves),

Foot switch (vacuum regulator),

Technical data:

Air flow rate up to 40 l/min,

The vacuum: - 0.9 bar at 686 mmHG,

Energy consumption 185 W,

Power supply 115/230 V, 50-60 Hz,

Dimensions 290 x 280 x 360 mm,

Net weight 8.5 kg,

He has a current Technical Passport issued,

Warranty 12 months

