



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NOUVAG Vacuson 60 suction pump / mammal (New)

New,
The Swiss production,
Advantages of the VACUSON 60 Mammal 60:
Economical thanks to excellent price/performance ratio,
Precise and infinitely adjustable suction capacity,
Robust, reliable and proven technology,
Low noise and vibration levels,
Dual piston pump without oil,
The Vacuson 60 suction pump motor drives two pistons to create a vacuum which allows an output of 60 litres per minute,
The Vacuson 60 is designed for suction and removal of fluids and intraoperative contents with adjustable suction force, which allows the device to be adjusted to the needs of the treatment,
The suction cup has a manometer to control the settings,
Piston technology used in Vacuson 60 mammals allows to obtain high vacuum parameters with pressure up to 0.9 bar and to work for many years without any decrease in parameters,
Thanks to their advantages, the devices can be used in both small and very large procedures,
The Vacuson 40 is equipped with an automatic overfilling system,
Standard equipment:
Control unit VACUSON 60 (pressure gauge 0-1 bar, vacuum regulator, 2 jar holders, 2 sockets: for ON/OFF footswitch, power cord, suction pipe with set of connection cables and disposable bacterial filter),
Footswitch On/Off,
Disposable bacterial filter (10pcs),
Power cable (5m),
Additional equipment (surcharge):
Jars of different capacities,
The "STAND" mammal trolley (2 shelves)
SANICAR mammal trolley (3 shelves),
Technical data:
Air flow rate up to 60 l/min,
The vacuum: - 0.9 bar at 686 mmHG,
Energy consumption 250 W,
Power supply 115/230 V, 50-60 Hz,
Dimensions 290 x 280 x 360 mm,
Net weight 12 kg,
He has a current Technical Passport issued,
Warranty 12 months

