



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DRAGER OXYLOG 1000 ambulance respirator - ventilator (Reconditioned)

Used, very good condition,
Dutch production,

The Oxylog 1000 is designed to ventilate emergency patients. Designed for use in the field, the Oxylog 1000 is one of the most versatile devices of its kind, with intuitive operation, robustness and mobility,

The ventilator has optical and acoustic alarms that control both airway and supply pressure, contributing to increased patient safety,

The power supply for all functions of the device is pneumatic, so the ventilator is independent of the power supply. Thanks to its design, the Oxylog 1000 is extremely reliable, which is so important to cope with the difficult conditions surrounding emergency ventilation,

Built-in manometer: -10 to 80 mbar,

PV indicator - supply pressure,

Adjustable ventilator functions:

Ventilation frequency (Freq): 5 - 35 1/min,

Volume One Minute (MV): 3 - 20 l/min,

Maximum airway pressure (Pmax): 25 - 55 mbar,

Mixing switch: 100% O2 or 60% O2,

Included:

Oxygen supply hose with a customer's choice of nozzle (AGA or DIN),

The patient's breathing hose,

PEEP valve valve,

Oxygen mask size 5,

User manual in Polish and English (PDF),

Dimensions: 22 x 21 x 12.5 cm,

Weight: 3.3 kg,

Cat 04,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan

