



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## ISO9001 Certified COVID-19 IgM/IgG Ab - Core tests (Cassette - Plate) (Family Package - Small) (New)

A new one,  
P.R.C. production,  
The rapid platelet test (Core COVID-19 IgM/IgG Ab Test) detects the presence and type of antibodies (IgM and IgG) produced by the immune system of a human coronavirus infected with SARS-CoV-2,  
No quick test is as effective as this. Its operation is supported by tests and certificates: medical CE and ISO 9001 and 13485. It has also been approved by the Polish Office for Registration of Medicinal Products (URPL).  
The COVID-19 IgM/IgG Ab is designed for professional in vitro diagnostics,  
The test is based on a detection method using colloidal gold, which allows for the determination of Covid-19 antibodies in IgM/IgG classes in patients,  
IgM antibodies can be detected by the test 4-12 days after infection. IgG antibodies are then detected,  
The packaging includes:  
COVID-19 IgM/IgG Ab test cassette - 5 pieces,  
Buffer - 5 pieces,  
A puppet - 5 pieces,  
Instructions for use in Polish - 1 piece,  
Storage temperature: 2 - 30°C,  
Set contains 5 pieces  
Specificity - 100%, false positive 0%.  
Sensitivity 97%, False negative 2.4%.  
No rapid test is as effective as this.  
The product has the highest effectiveness, confirmed by tests and certificates  
The last picture shows the original box.

