



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ISO9001 certified COVID-19 IgM/IgG Ab - Core tests (cassette - plate) (1 pc) (New)

A new one,
P.R.C. production,
The rapid platelet test (Core COVID-19 IgM/IgG Ab Test) detects the presence and type of antibodies (IgM and IgG) produced by the immune system of a human coronavirus infected with SARS-CoV-2,
No quick test is as effective as this. Its operation is supported by tests and certificates: medical CE and ISO 9001 and 13485. It was also approved by the Polish Office for Registration of Medicinal Products (URPL),
The COVID-19 IgM/IgG Ab is designed for professional in vitro diagnostics,
The test is based on a detection method using colloidal gold, which allows for the determination of Covid-19 antibodies in IgM/IgG classes in patients,
IgM antibodies can be detected by the test 4-12 days after infection. IgG antibodies are then detected,
The package includes:
COVID-19 IgM/IgG Ab test cassette - 1 piece,
Buffer - 1 piece,
A puppet - 1 piece,
Instructions for use in Polish - 1 piece,
Storage temperature: 2 - 30°C,
Set contains 1 piece
Specificity - 100%, false positive 0%.
Sensitivity 97%, False negative 2.4%.
No rapid test is as effective as this.
The product has the highest effectiveness, confirmed by tests and certificates
The last picture shows the original box.

