



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **GE STENOSCOPI C arm X-ray device (Reconditioned)**

Used very good condition,  
USA production,  
Date of manufacture: 2000 r.  
Power supply: 240 V,  
Frequency: 50/60 Hz,  
Very easy to use, pressing the Stenoscopy foot pedal displays the perfect photo, automatically adjusting to the overexposed part of the body,  
Specially designed shape and collimator reduce external radiation,  
The function measuring the amount of radiation received by the skin allows to control the radiation dose,  
Electric arm height adjustment,  
Control panel equipped with:  
Two 19" screen diagonal monitors,  
Two handles to move,  
Mobile on wheels with brakes,  
SONY Model UP980CE printer:  
Used, very good condition,  
Japanese production,  
Thermal printing with 1280 points and 256 grey shades,  
Multifunctional printing,  
Print speed: CCIR: approx. 12 seconds EIA: approx. 10 seconds,  
Accepts both 15.5 kHz and 31.5 kHz composite video signal,  
SPLIT function - printing 1, 2, 4 or 6 images on one page,  
The brightness and contrast can be adjusted,  
Designed for printing CCTV systems and medical equipment such as ultrasound, X-ray,  
Auto cut - cuts paper automatically when printing is completed,  
It also prints images on blue transparent film, Sony UPT-210BL,  
Paper types: Sony UPP-210HD or Sony UPP-210SE,  
The camera is after inspection, with important acceptance tests,  
Current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (Poland only): Instalments, Leasing, Loan  
OVERVIEW PICTURES OF THE DEVICE

