



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Drip stand (white) (Reconditioned)

Used, very good condition,
German production,
Tripod made of powder-coated steel,
The mast is equipped with 2 hangers,
Infinitely variable height adjustment in range: 125 - 218 cm,
The round base has 3 legs,
Base diameter: 39 cm,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

