



## **WEINMANN MEDUMAT VARIABEL respirator - ventilator (Cat. 03) (Reconditioned)**

Used, very good condition,  
German production,  
Infinitely variable flow parameters,  
Built-in adjustable flow meter module,  
Equipped with a manometer,  
The ventilator has a mammal function,  
Built-in inhaler module,  
The ventilator is placed in a special metal suitcase,  
Included:  
Korgiel W8 O2 oxygen reducer,  
Oxygen supply hose for AGA-type cylinders,  
Patient breathing hose with valve and PEEP valve,  
Patient mask size 5,  
Nasal cannula with a 2.1 m drain,  
Medical oxygen cylinder with capacity: 2 l,  
Dimensions: 48.5 x 15 x 43 cm,  
Weight: 12.6 kg,  
Cat 03,  
He has a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (Poland only): Instalments, Leasing, Loan

