



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



EUROMI fat transfer kit NLF1 (Fat Collector / Fat Transfer) (New)

A new one,
Belgian production,
NLF1 allows the separation of fat from serum by decantation,
The set includes:
Sterile 1500 ml fat container,
Single-use pass-through sampling cannula consisting of two parts (VITERBO cannula, diameter 3 mm, length of the working part of the cannula - 250 mm, length of the whole cannula - 540 mm),
Suction hose with buckle,
Pressure hose,
Hose with clamp for taking fat from the container,
Two elbows (attached to the air outlets of the head) to prevent the moving cannula from rubbing against adjacent wires (photo 8),
Reusable stand made of stainless steel,
The container is made of high-quality, unbreakable and impact-resistant plastic,
It has clear liquid capacity markings,
The system is biocompatible,
Cytotoxicity certificates,
The user manual in English,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan,
He has a current Technical Passport issued

