



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



GE VIVID USG 7 (Reconditioned)

Used, good technical condition,
U.S.A. production,
High class clinical apparatus supporting electronic broadband probes: sector, transesophageal, intraoperative, convex, microconvex, linear and
2D and 3D/4D matrix probes,
Imaging modes:
2D,
M-Mode,
High Frame Rate Color Flow,
Color Doppler Velocity,
Color Flow Mapping,
Color Doppler Energy,
Color Power Doppler,
Directional Tissue Imaging (DTI),
Pulse Wave,
Continuous Wave (CW),
3-dimensional multi-dimensional imaging, spatial rendering,
Stress Echo tests with quantitative analysis with the possibility to include contrast tests,
The device is controlled in an easy and intuitive way,
Programmable functions and reports,
Functions:
B-FLOW,
SRI,
DICOM,
Possibility to write images on CDs in formats: BMP, TIFF, AVI, MPEG,
Equipped with a 21" color LCD monitor with tilt and swivel capability,
The extremely ergonomic design of the device, together with the possibility of height and angle adjustment of the console, makes it very easy
to adjust it to the user's needs,
Central brake,
Included:
Linear head GE 10L,
GE sector head M3S,
Videoprinter SONY UP-895MDW,
Power cord,
User manual in English (PDF),
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

