



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SIEMENS SERVO VENTILATOR 300 respirator - ventilator (Reconditioned)

Used, very good condition,
Swedish production,
Number of hours worked: 58003 h,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 110 W,
Built-in digital displays,
Working modes: adults, children, newborns,
Adjustable ventilator functions:
CMV freq: 10 - 150 bpm,
Insp time: 10 - 80%,
Pause time: 0 - 30%,
Insp. Rise time: 0 - 10%,
SIMV freq: 0 - 40 b/min,
O2 conc.: 21 - 100%,
Press. Control level above PEEP: 0 - 100,
Press. Support level above PEEP: 0 - 100,
PEEP: 0 - 50,
Trig. Sensitivity Level below PEEP: -16 - 0,
It has acoustic alarms,
Foiled control panel,
Diode bar displays,
Rails for fixing additional accessories,
Included:
Patient's respiratory system with a 2 l breathing bag,
Oxygen (O2) and air (AIR) supply hoses with AGA-type connectors,
Power cord,
User manual in English (PDF),
Respirator placed on a mobile trolley on four large wheels, including two with brakes,
Overall dimensions:
Width: 54 cm,
Depth: 70 cm,
Height: 118.5 cm,
Weight: 43 kg,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

